



A & E Lab (UK) Co., Ltd

Hotplates

Application:

We offer a variety of hotplates, ranging from different size, different top materials, and different temperature setting. Most hotplates come with the spill prevention feature, and hot warning features.

Technical Specifications:

Model	AHC-101 Microcrystal Top Hotplates	AHC-102 Microcrystal Top Hotplates
Plate Dimension	160X160mm	300X300mm
Heated Area	120X120mm	220X220mm
Heated Power	550W	1500W
Maximum Plate Temp.	550°C	550°C
Dimension (mm,)	220X240X145	330X380X145
Net Weight (Kg)	2.8	8
Power Supply	230V, 50Hz or 110V, 60Hz	



AHC-101/102

Model	AHA-101 Aluminum Top Hotplates	AHA-102 Aluminum Top Hotplates
Plate Dimension	127(5')x127(5')mm	150 Diameter Round Top
Heated Area	100X100mm	120 Diameter
Heated Power	350W	350W
Maximum Plate Temp.	350°C	350°C
Dimension (mm)	127X130X85	190X260X125
Net Weight (Kg)	1.0	2.0
Power Supply	230V, 50Hz or 110V, 60Hz	



AHA-101



AHA-102

