

CAPP

D E N M A R K

 **ecopipette**

SINGLECHANNEL PIPETTES

ecopipette.capp.dk


ecopipette Innovation

- THIS PIPETTE IS BUILT TO LAST. MAKING NO COMPROMISE ON QUALITY IN PRODUCT DEVELOPMENT, RAW MATERIALS USED AND MANUFACTURING PROCESSES WE DEVELOPED THE GREENEST PIPETTE AVAILABLE ON THE MARKET. THE **ECOPIPETTE** IS NOT ONLY PACKAGED USING RECYCLABLE, BIODEGRADABLE MATERIALS, BUT CONSTRUCTED FROM SOME OF THE MOST RENEWABLE DURABLE RESOURCES AVAILABLE TO ENSURE AN ADVANCED LIFE SPAN AND GREENEST AFTER SERVICE RECYCLING.
- CAPP ECOTRADE™ PROGRAM: RETURN ANY BRAND OF OLD PIPETTES IN EXCHANGE FOR DISCOUNTED NEW ECOPIPETTES™. AT CAPP APS, WE WILL ENSURE AN ENVIRONMENTALLY SAFE WASTE DISCHARGE OF ALL OLD PIPETTES. ASK YOUR FAVOURITE AUTHORIZED CAPP DISTRIBUTOR FOR MORE INFORMATION ABOUT OUR ECOTRADE™ PROGRAMME.

CAPP **ECOPIPETTE**™ UNIQUE ADVANTAGES:

- LOCKABLE VOLUME CONTROLLER REQUIRING 50% FEWER TURNS
- DISPLAY IS FULLY VISIBLE DURING HAND OPERATIONS
- REDUCED FORCE REQUIREMENT FOR ASPIRATION AND BLOW-OUT
- THE ONLY PIPETTE WITH FIXED AND VARIABLE EXCHANGEABLE COUNTER MECHANISM
- LARGE PAD FOR EASE OF TIP EJECTION
- FULLY AUTOCLAVABLE PIPETTE BODY (NO DISASSEMBLY REQUIRED)
- DOUBLE O-RING FOR OPTIMAL TIP SEAL
- COMPATIBLE WITH ALL MAJOR PIPETTE TIP BRANDS
- ROBUST DESIGN INCLUDES METAL TIP BASE
- LOWEST INCIDENCE OF BREAKAGE RESULTING FROM ACCIDENTAL DROPPING OR TIPS MOUNTING.



 CAPP APS
 LANDBRUGSVEJ 10
 DK-5260 ODENSE
 DENMARK

TEL. (+45) 6613 6140
 FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK
 WEB: WWW.CAPP.DK



Ordering Info

CAT. NO.	DESCRIPTION	INACCURACY %	IMPRECISION %	COLOR
C02-1S	CAPP PIPETTE, VARIABLE VOL. 0.2-2 UL, ECOPIPETTE	12.00 /2.00	6.00 /1.00	
C10-1	CAPP PIPETTE, VARIABLE VOL. 0.5-10 UL, ECOPIPETTE	2.50 /1.00	1.80 /0.50	
C20-1	CAPP PIPETTE, VARIABLE VOL. 2-20 UL, ECOPIPETTE	5.00 /1.00	1.50 /0.30	
C50-1	CAPP PIPETTE, VARIABLE VOL. 5-50 UL, ECOPIPETTE	2.00 /1.00	1.20 /0.30	
C100-1	CAPP PIPETTE, VARIABLE VOL. 10-100 UL, ECOPIPETTE	2.00 /0.70	0.80 /0.20	
C200-1	CAPP PIPETTE, VARIABLE VOL. 20-200 UL, ECOPIPETTE	1.50 /0.60	0.50 /0.15	
C1000-1	CAPP PIPETTE, VARIABLE VOL. 100-1000 UL, ECOPIPETTE	1.50 /0.80	0.40 /0.15	
C5000-1	CAPP PIPETTE, VARIABLE VOL. 1000-5000 UL, ECOPIPETTE	1.20 /0.60	0.30 /0.15	

STARTER PACK

SKE-03-XX THREE PIPETTES FREE OF CHOICE, STAND, RACKS OF TIPS AND TIMER